

FAPAS QC MATERIAL DATA SHEET	T3106QC
Matrix	Beverage
Weight / Volume of Contents	2 x 200 ml

Analyte	Qualitative Results
Alcohol (Ethanol) TMA	POSITIVE (70% consensus) Advisory level 0.01% alcohol (ethanol) by volume
Alcohol (Ethanol) TMB	POSITIVE (91% consensus) Advisory level 0.05% alcohol (ethanol) by volume

This data sheet is applicable until	02 Jul 2018
Recommended Storage on receipt	-20°C

Notes

- Mix the QC material thoroughly before taking a representative analytical sample.
- Stability of the QC material has been established as sufficient for the scope of the
 proficiency test from previous experience, expert advice and published literature. FAPAS
 advises that the QC material is analysed within the recommended date. FAPAS QC
 materials are intended to be used as single-analysis samples.
- Full details on the proficiency test procedure used to characterise this QC material are available in the Protocol, Part 1 - Common Principles, freely available to download from the FAPAS website.