



Fera Science Ltd. (Fera)
 Sand Hutton, York, UK, YO41 1LZ
 Tel: +44 (0)1904 462100 Fax: +44 (0)1904 500440
 info@fapas.com www.fapas.com

Fapas® QC MATERIAL DATA SHEET	T02376bQC
Matrix	Blank Bovine Urine
Weight / Volume of Contents	25 ml

Analyte	Assigned Value, x_a	Range for $ z \leq 2$	Units	No. of data points producing x_a
Blank Material	No Assigned value set			

This data sheet is applicable until	24 Feb 2021
Recommended Storage on receipt	-20°C

Notes

- This is a blank test material and cannot be purchased on its own; this test material should only be used in conjunction with T02376QC veterinary drug residue material.
- Mix the QC material thoroughly before taking a representative analytical sample.
- This blank test material has been screened to a level of 0.6 µg/kg, this does not preclude the ability of sensitive methods to find pesticide residues below this level.
- Fapas® advises that the QC material is analysed within the recommended date. Fapas® QC materials are intended to be used as single-analysis samples.
- Full details on the proficiency test procedure used to characterise this QC material are available in the Protocol, Part 1 - Common Principles, freely available to download from the Fapas® website.