



Fapas® QC MATERIAL DATA SHEET	T02518bQC
Matrix	Blank Chicken (Hens) Eggs
Weight / Volume of Contents	20 g

Analyte	Assigned Value, x_a	Range for $ z \leq 2$	Units	No. of data points producing x_a
Blank for analytes in T02518QC	No Assigned value set			

This data sheet is applicable until	31 Aug 2025
Recommended Storage on receipt	-20°C

Notes
<ul style="list-style-type: none">• This is a blank test material and cannot be purchased on its own; this test material should only be used in conjunction with T02518QC veterinary drug residue material.• Mix the QC material thoroughly before taking a representative analytical sample.• This blank material was screened for the presence of any incurred coccidiostats. Any very low level detections will not significantly affect the use of this blank material for matrix-matched calibrations or recoveries at the concentrations of the associated QC sample.• Fapas® advises that the QC material is analysed within the recommended date. Fapas® QC materials are intended to be used as single-analysis samples.• Full details on the proficiency test procedure used to characterise this QC material are available in the Protocol, Part 1 - Common Principles, freely available to download from the Fapas® website.