



Fapas® QC MATERIAL DATA SHEET	T05174bQC
Matrix	Blank Olive Oil (Virgin)
Weight / Volume of Contents	50g

Analyte	Assigned Value, x_a	Range for $ z \leq 2$	Units	No. of data points producing x_a
Blank for analytes in T05174QC	No Assigned value set			

This data sheet is applicable until	26 Oct 2026
Recommended Storage on receipt	-20°C

Notes
<ul style="list-style-type: none">• This is a blank test material and cannot be purchased on its own; this test material should only be used in conjunction with T05174QC pesticide residue material.• Mix the QC material thoroughly before taking a representative analytical sample.• This blank test material has been screened to a level of 20 µg/kg, this does not preclude the ability of sensitive methods to find pesticide residues below this level.• Fapas® advises that the QC material is analysed within the recommended date. Fapas® QC materials are intended to be used as single-analysis samples.• Full details on the proficiency test procedure used to characterise this QC material are available in the Protocol, Part 1 - Common Principles, freely available to download from the Fapas® website.• The levels of residues present should be determined in the form it is received (no further correction for wet weight or dry weight or reconstitution).