



|                               |             |
|-------------------------------|-------------|
| Fapas® QC MATERIAL DATA SHEET | T02581bQC   |
| Matrix                        | Blank Honey |
| Weight / Volume of Contents   | 40 g        |

| Analyte                        | Assigned Value, $x_a$ | Range for $ z  \leq 2$ | Units | No. of data points producing $x_a$ |
|--------------------------------|-----------------------|------------------------|-------|------------------------------------|
| Blank for analytes in T02581QC | No Assigned value set |                        |       |                                    |

|                                     |             |
|-------------------------------------|-------------|
| This data sheet is applicable until | 20 Dec 2026 |
| Recommended Storage on receipt      | -20°C       |

|   |
|---|
| Notes   |
| <ul style="list-style-type: none"><li>• This is a blank test material and cannot be purchased on its own; this test material should only be used in conjunction with T02581QC veterinary drug residue material.</li><li>• Mix the QC material thoroughly before taking a representative analytical sample.</li><li>• Fapas® advises that the QC material is analysed within the recommended date. Fapas® QC materials are intended to be used as single-analysis samples.</li><li>• Full details on the proficiency test procedure used to characterise this QC material are available in the Protocol, Part 1 - Common Principles, freely available to download from the Fapas® website.</li><li>• The test material was screened for the presence of any incurred analytes. None were found at or above the normal reporting limit of the sub-contracted analytical laboratory.</li></ul> |