



Fapas® QC MATERIAL DATA SHEET	T02585bQC
Matrix	Blank Honey
Weight / Volume of Contents	20 g

Analyte	Assigned Value, $x_a$	Range for $ z  \leq 2$	Units	No. of data points producing $x_a$
Blank for analytes in T02585QC	No Assigned value set			

This data sheet is applicable until	09 Mar 2027
Recommended Storage on receipt	-20°C

Notes
<ul style="list-style-type: none"><li>• This is a blank test material and cannot be purchased on its own; this test material should only be used in conjunction with T02585QC veterinary drug residue material.</li><li>• Mix the QC material thoroughly before taking a representative analytical sample.</li><li>• Fapas® advises that the QC material is analysed within the recommended date. Fapas® QC materials are intended to be used as single-analysis samples.</li><li>• Full details on the proficiency test procedure used to characterise this QC material are available in the Protocol, Part 1 - Common Principles, freely available to download from the Fapas® website.</li><li>• The test material was screened for the presence of any incurred analytes. None were found at or above the normal reporting limit of the sub-contracted analytical laboratory.</li></ul>