



Fapas® QC MATERIAL DATA SHEET	T3504AQC
Matrix	Flour
Weight / Volume of Contents	125 g

Analyte	Qualitative Results
Insect	NOT DETECTED (86% consensus)
Metal	NOT DETECTED (100% consensus)
Glass	NOT DETECTED (100% consensus)
Plastic	DETECTED (100% consensus)
Stones	DETECTED (71% consensus)
Wood	NOT DETECTED (100% consensus)
Other	NOT DETECTED (83% consensus)

This data sheet is applicable until	20 Mar 2027
Recommended Storage on receipt	Ambient

Notes
<ul style="list-style-type: none"><li>• Mix the QC material thoroughly before taking a representative analytical sample.</li><li>• Stability of the QC material has been established as sufficient for the scope of the proficiency test from previous experience, expert advice and published literature. Fapas® advises that the QC material is analysed within the recommended date. Fapas® QC materials are intended to be used as single-analysis samples.</li><li>• Full details on the proficiency test procedure used to characterise this QC material are available in the Protocol, Part 1 - Common Principles, freely available to download from the Fapas® website.</li></ul>