



| | |
|-------------------------------|----------------------|
| Fapas® QC MATERIAL DATA SHEET | T19445bQC |
| Matrix | Blank Cucumber Purée |
| Weight / Volume of Contents | 90 g |

| Analyte | Assigned Value, x_a | Range for $ z \leq 2$ | Units | No. of data points producing x_a |
|--------------------------------|-----------------------|------------------------|-------|------------------------------------|
| Blank for analytes in T19445QC | No Assigned value set | | | |

| | |
|-------------------------------------|-------------|
| This data sheet is applicable until | 10 Mar 2027 |
| Recommended Storage on receipt | -20°C |

| |
|---|
| Notes |
| <ul style="list-style-type: none">• This is a blank test material and cannot be purchased on its own; this test material should only be used in conjunction with T19445QC pesticide residue material.• Mix the QC material thoroughly before taking a representative analytical sample.• Fapas® advises that the QC material is analysed within the recommended date. Fapas® QC materials are intended to be used as single-analysis samples.• Full details on the proficiency test procedure used to characterise this QC material are available in the Protocol, Part 1 - Common Principles, freely available to download from the Fapas® website.• The test material was screened for the presence of any incurred analytes. None were found at or above the normal reporting limit of the sub-contracted analytical laboratory. |