



Fapas® QC MATERIAL DATA SHEET	T19452bQC
Matrix	Blank Coffee (Green, Milled)
Weight / Volume of Contents	90 g

Analyte	Assigned Value, x_a	Range for $ z \leq 2$	Units	No. of data points producing x_a
Blank for analytes in T19452QC	No Assigned value set			

This data sheet is applicable until	08 Feb 2027
Recommended Storage on receipt	-20°C

Notes
<ul style="list-style-type: none">• This is a blank test material and cannot be purchased on its own; this test material should only be used in conjunction with T19452QC pesticide residue material.• Mix the QC material thoroughly before taking a representative analytical sample.• Fapas® advises that the QC material is analysed within the recommended date. Fapas® QC materials are intended to be used as single-analysis samples.• Full details on the proficiency test procedure used to characterise this QC material are available in the Protocol, Part 1 - Common Principles, freely available to download from the Fapas® website.• The test material was screened for the presence of any incurred analytes. None were found at or above the normal reporting limit of the sub-contracted analytical laboratory.