



FAPAS QC MATERIAL DATA SHEET	T05123bQC
Matrix	Blank Olive Oil
Weight / Volume of Contents	50 g

Analyte	Assigned Value, $x_a$	Range for $ z  \leq 2$	Units	No. of data points producing $x_a$
Blank Material	No Assigned value set			

This data sheet is applicable until	02 Nov 2020
Recommended Storage on receipt	Ambient

Notes
<ul style="list-style-type: none"><li>• This is a blank test material and cannot be purchased on its own; this test material should only be used in conjunction with T05123QC pesticide residue material.</li><li>• Mix the QC material thoroughly before taking a representative analytical sample.</li><li>• This blank test material has been screened to a level of 10 µg/kg, this does not preclude the ability of sensitive methods to find pesticide residues below this level.</li><li>• FAPAS advises that the QC material is analysed within the recommended date. FAPAS QC materials are intended to be used as single-analysis samples.</li><li>• Full details on the proficiency test procedure used to characterise this QC material are available in the Protocol, Part 1 - Common Principles, freely available to download from the FAPAS website.</li></ul>